Leica FS M

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Holger Grasse

(Product Compliance Manager R&D)

In case of questions, please contact the hotline:

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Leica
MICROSYSTEMS

Copyrights

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Nevertheless, no liability of any kind may be assumed for the accuracy of this manual's contents. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information in this manual is subject to modification at any time and without prior notification.

Contents

1. Important Notes about this Manual.....	6	9. Accessories	21
2. Intended Use of the Device	7	9.1 Motorized work table.....	21
3. Safety Notes.....	8	9.2 Stackable stages.....	22
3.1 General safety notes.....	8	9.3 Slide holder.....	24
3.2 Electrical safety.....	8	9.4 Universal rotating mechanism	31
3.3 Disposal.....	9	9.5 Light panels	31
4. Installation Site.....	10	9.6 Filter for color contrast.....	32
5. Unpacking	11	9.7 Further illumination options.....	33
6. Overview of the Instrument.....	13	9.8 Stage scales and calibration standard.....	34
7. Assembly	14	10. Technical Data.....	35
7.1 Comparison bridge.....	14	11. Care of the Instrument.....	38
7.2 Illumination	14	11.1 Dust cover.....	38
7.3 Objectives.....	15	11.2 Cleaning	38
7.4 Tube	15	11.3 Handling acids and alkalines	39
7.5 Eyepieces.....	15	12. Essential Wear and Spare Parts	40
7.6 Connecting to the mains	16	13. EC Declaration of Conformity	41
8. Operation	17		
8.1 Switching on.....	17		
8.2 Stages.....	17		
8.3 Tube	17		
8.4 Eyepieces.....	18		
8.5 Brightness adjustment.....	18		
8.6 Focusing.....	18		
8.7 Objectives.....	18		
8.8 Magnification changer	19		
8.9 Adjusting the height of the stand.....	19		
8.10 Diaphragm setting	19		
8.11 Magnification calibration.....	20		

1. Important Notes about this Manual

Caution!

This operating manual is an essential component of the comparison microscope Leica FSM and must be read carefully before the microscope is put into operation or used.

This operating manual contains important instructions and information for the operational safety and maintenance of the comparison microscope and accessories. Therefore, it must be kept and taken care of.

Text symbols and their meanings:

(1.2)

Numbers in parentheses such as (1.2), correspond to illustrations (in the example, Figure 1, Item 2).

→ p.20

Numbers with pointer arrows (for example → page 20), point to a certain page of this manual.

Special safety instructions are indicated with the triangle symbol shown here, and have a gray background.

Caution! The microscope and accessories can be damaged when operated incorrectly.

Danger due to hot surface.

Instructions for disposing of the instrument, its accessories and consumables.

Explanatory note.

*

Item not contained in all configurations.

2. Intended Use of the Device

Comparisons of small surfaces or profiles need to be performed with sufficient precision, if the specimens of comparison can be observed simultaneously.

If only a macroscope with just one imaging beam path is available for such comparisons, then, depending on the interruption of observation by a change of objective, high demands are made on the memory of the observer. Therefore, it is impossible to rule out incorrect interpretation in the case of specimens with very minor differences in the surface or profile.

The Leica FS M Comparison Macroscope with two imaging beam paths, is designed for comparing structures of items such as those listed below in split images and composite images in forensic laboratories, document printing facilities, mints and banks, and documenting the results on photographs and videos:

- Evidence on fired ammunition pieces
- Toolmarks
- Documents

Many types of illuminators are available for different applications (including UV light).

The instrument is compliant with EC Directives 2006/95/EG on low voltage electrical equipment and 2004/108/EG regarding electromagnetic compatibility for laboratory use.

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from using the device for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In the event of such, the declaration of conformity shall become null and void.

3. Safety Notes

3.1 General safety notes

This safety class 1 instrument was built and tested in accordance with:

EN 61010-1:2001,

IEC 61010-1:2001,

Safety requirements for electrical equipment for measurement, control and laboratory use.

Caution!

To maintain the unit in its original condition and to ensure safe operation, the user must follow the instructions and warnings contained in this User Manual.

Caution!

The instruments and accessories described in this operating manual have been tested for safety and potential hazards.

The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the instrument is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the instrument or noncompliant use shall void all rights to any warranty claims and product liability

3.2 Electrical safety

General specifications

For indoor use only.

Supply voltage: 100-240 VAC

Frequency: 50-60 Hz

Power input: I_{max} 0.5 A

Fuse (power supply): 2 x 2,5 A H 250 V

Ambient temperature: 15-35°C

Relative humidity: max. 80% up to 30°C

Overvoltage category: II

Pollution degree: 2

Caution!

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the instrument, or release of the ground wire connection, can cause the instrument to become hazardous. Intentional ground interruption is not permitted!

Type label Leica FS M

Caution!

By connecting them to ground, ancillary devices providing their own and/or separate power supplies and which are connected to the comparison microscope can be connected to the same ground potential. For connections without a ground connector, Leica Service must be consulted.

Caution!

Never use any fuses as replacements other than those of the types and the current ratings listed here. Bypassing the fuse holder is not permitted.

Caution!

The comparison microscope's electrical accessory components are not protected against water. Water can cause electric shock.

Do not assemble the comparison microscope and its accessory components in the immediate vicinity of a water connection or other areas where the possibility of water ingress exists.

Caution!

Protect the comparison microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.

Caution!

Before exchanging the fuses or lamps, be absolutely certain to switch-off the main power switch and remove the power cable.

3.3 Disposal

Once the product has reached the end of its service life, please contact Leica Service or Sales about disposal.

Please observe and comply with the national and federal laws and regulations that are equivalent to EC directives such as WEEE.

Note!

Like all electronic devices, this instrument, its accessory components and consumables must never be disposed of with general household waste.

4. Installation Site

Caution!

Ensure that the workplace area is free of oil and chemical fumes. Vibrations, direct sunlight and strong temperature fluctuations interfere with measurements as well as microphotography and image documentation.

The basic prerequisite is a stable instrument table of optimum height (70–80 cm). Combined with an ergonomically correct, adjustable chair, these are the basic requirements for fatigue-free microscopy.

Allowable ambient conditions

Temperature	15–35°C
Relative humidity	max. 80% up to 30°C

Macrosopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

For additional instructions, see the chapter entitled "Care of the Instrument" → p. 38.

Caution!

Electrical components must be assembled at least 10 cm from the wall and away from flammable substances.

Caution!

When installing the macroscope, make sure the mains socket is freely accessible so that the instrument can be quickly disconnected from the mains if necessary.

5. Unpacking

The stand is delivered in a solid box with packaging inserts.

The delivery includes the following components:

- Stand with base, as well as the stages, which are already installed. The stand is equipped with two carrying handles.
- Comparison bridge
- Tube
- One or two cold light sources
- Accessory case with slide holders
- Optical accessories
- Oblique incident illuminator
- Optional accessories

Fig. 1 Stand with carrying handles
1 Removable carrying handle

When unpacking the equipment, follow these instructions:

- After opening the box, first remove the comparison microscope and then the cardboard boxes, which hold the accessories.
- Lift the stand out of the carton by its carrying handles (1.1).

Caution!

When lifting the instrument out of the packaging, and when transporting it to its installation site, carry it by its handles only.

Caution!

Note that the stand weighs approx. 32 kg. Two people are required to lift it out of its packaging and carry it to its installation site.

Later, the carrying handles can be unscrewed and removed from the upper part of the microscope. Then, seal the threaded holes using seal plugs.

5. Unpacking

- Carefully remove all components from the transportation and packaging materials.

Caution!

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces of optical surfaces in a short time. For additional instructions, see the chapter entitled "Care of the Instrument" → p. 38.

- Check the delivery carefully against the packing list, delivery note or the invoice. We highly recommend keeping a copy of these documents with the manual in order to have delivery date and delivery package information available when (for example) placing repeat orders or making service calls. Please make sure that no small parts are lost in the packaging materials.

Note:

Keep the packaging materials for storage purposes and transportation of the individual components.

Caution!

Do not yet connect the comparison microscope and peripherals to the power supply at this point!

6. Overview of the Instrument


Fig. 2 Leica FS M

- | | |
|---|--|
| <ul style="list-style-type: none"> 1 Comparison bridge 2 Leica DFC digital camera 3 Tube with eyepieces 4 Knob for setting bridge mode, right 5 Positioning of dividing line 6 Adjustment screws for dividing line 7 Objective turret with objectives 8 LED spot illumination | <ul style="list-style-type: none"> 9 Coarse and fine focusing, right 10 Synchronized movement of both stages in x direction 11 Rocker switch for motorized adjustment of comparison bridge height 12 Coarse and fine focusing, left 13 Specimen stage with coaxial drive 14 Mount for filter slides 15 Knob for setting bridge mode, left |
|---|--|

7. Assembly

Only a few commonly used screwdrivers and keys are necessary for assembly; these are included in the delivery package.


There is a holder on the back of the stand for keeping the dividing line adjustment key.

! Caution!

Unscrew the transit protection screws (3.1) on the underneath of the comparison bridge and screw them into the threads provided for safe-keeping (3.2) (parking position).

Fig. 3 Transit protection for the comparison bridge

- 1 Transit protection screws
- 2 Parking position for the transit protection screws

7.1 Comparison bridge

- From above, set the comparison bridge onto the dovetail guide and push the bridge downwards until the upper end of the dovetail guide is flush.
- Fasten the comparison bridge using the side clamping screw.

7.2 Illumination

If available:

- To accommodate the illuminator, assemble the rotating mechanism with articulated arm (Fig. 4) to the stage.

See Chapter "9. Accessories" → p. 31, 33 for further illumination options.

Fig. 4 Illumination unit

- 1 LED spot lamp
- 2 Articulated holder
- 3 Rotating mechanism

7.3 Objectives

The position of the individual objectives is preset at the factory. The object turret eyes are labeled (5.2).

! Caution!

The objectives for the left and right objective nosepiece are packed separately. When screwing in the objectives, make sure you put them in the right side.

- Screw the objectives into the objective turret.

Position	Objective
50	0.4x
100	1x
200	2x
500	4x
1000	8x

Note:

The SL 0.3:1 auxiliary lens (order No. 11 532 301), which is attached to the objective, is available for the objective 0.4x. 1:1 imaging is achieved using the auxiliary lens.

7.4 Tube

- Slacken the clamp screw on the tube mount
- Put the tube into the mount.
- Retighten the clamp screw.

Note:

If using the optional magnification changer, insert it between the tube mount and the tube and fix in position with the clamp screw.

7.5 Eyepieces

- The eyepieces are inserted into the eyepiece tubes on the tube.

Fig. 5 Objective turret with objectives

- 1 Filter receptacle
- 2 Objective positions

7. Assembly

7.6 Connecting to the mains

- Plug the power cord into the back of the comparison microscope (6.1).
- Then connect the comparison microscope to the mains.

Fig. 6 Back view of the Leica FS M

1 Mains connection

8. Operation

The Leica FS M Comparison Microscope offers the following functions:

- Simultaneous (in x direction) or separate methods for the left and right stage
- Left and right focusing
- Left and right objective nosepiece
- Adjustable stand height

8.1 Switching on

The stand is switched on with the switch (7.1). The illumination is switched on separately.

8.2 Stages

The stages on the left and right side can be moved either independently of each other with the left or right coaxial drive (7.2) or simultaneously in the x direction (7.3).

8.3 Tube

Note:

Close any unused tube openings, as otherwise stray light can interfere with observation.

Adjusting the viewing distance

- Adjust the viewing distance of the eyepiece tubes so that a congruent total image is seen.

For the Ergotube:

Adjusting the Angle of Vision

- The viewing angle can be adjusted from 5° to 32° by tilting the binocular eyepiece.

Fig. 7 Right side of the microscope

- 1 On/off switch
- 2 Coaxial drive
- 3 Synchronized movement of both stages in x direction

8. Operation

8.4 Eyepieces

Note:

The eyepiece's aperture protector must be removed during microscopy while wearing eye-glasses, or at least folded back.

Eyeglasses with multifocal lenses (bifocals and smooth view glasses) must be removed while operating the microscope.

Correction for Vision Problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.

If both oculars have an adjustable eyelens, the eyelens of the right ocular can be brought to the position for normal vision, which is indicated by an encircling line.

When using a graticule in the right eyepiece, the eyelens must first be focused on the graticule.

- Then, view the same specimen using the left eye and rotate the eyelens of the left eyepiece until the specimen appears in sharp focus. While doing so, do not use the focus dial.

8.5 Brightness adjustment

To set the intensity of the illumination, refer to the separate "LED Controller Coolcontrol" manual provided.

Note:

The Leica FS M Comparison Microscope offers a wealth of illumination options ranging from LED illumination, cold light illumination with fiber optic guides to light panels for diffuse illumination and special illumination units for UV light.

See also "9. Accessories" → p. 31, 33 and the corresponding special instructions.

8.6 Focusing

The coarse and fine focusing is operated for each side separately with the left and right focus knobs (8.1 and 8.3).

8.7 Objectives

The desired objectives are swiveled into place manually.

All objectives have an iris diaphragm. With the iris diaphragm, the NUM aperture can be adjusted continuously. Five ratchets make it easy to adjust to identical settings on both sides.

Fig. 8 Front view of the Leica FS M

- 1 Coarse and fine focusing, left
- 2 Rocker switch for setting the stand height
- 3 Coarse and fine focusing, right

8.8 Magnification changer*

The magnification changer is switched on and off by hand.

Magnification factors: 1x, 1.5x, 2x

When the magnification changer is switched on, the current magnification of both images is increased by the selected factor.

The magnification changer controls the visual and photo outputs simultaneously.

8.9 Adjusting the height of the stand

The rocker switch (8.2) adjusts the height of the stand individually. This is provided for aligning the height of specimens or specimen holders which vary in height.

8.10 Diaphragm setting


The diaphragms allow transition from a superimposed image to the separate view of the right and left image, as well as the setting of the dividing line.

Four diaphragm positions are already predefined and can be selected via the turning knobs (9.1 and 9.2):

- Mix** The left and right images are superimposed (composite image method). Both knobs at position ○
- LR** The left and right image appear side-by-side. A line in the middle of the field of vision separates the two images (split-image method). Left knob at position ●
Right knob at position ●
- L** Only the left image is displayed. Left knob at position ●
Right knob at position ●
- R** Only the right image is displayed. Left knob at position ●
Right knob at position ○

Fig. 9 Bridge mode controls

- 1 Diaphragm setting, left
- 2 Diaphragm setting, right
- 3 Screw for adjusting the dividing line (parallelism)
- 4 Positioning of dividing line
- 5 Center position for accessing the adjustment screws
- 6 Screw for adjusting the dividing line (width)

8. Operation

The dividing line is moved to the left or to the right by rotating the knob (9.4) and has an impact on the size of the fields in split image mode. The center position (9.5) is used to bring the dividing line back into the center.

The width of the dividing line is changed with the adjustment screw (9.6). The preferential setting should be a thin, black hairline to highlight the distinction between the two fields. If the edges of the dividing line do not run parallel to each other, parallelism can be set with the adjustment screw (9.3).

To access the adjustment screws, the pointer on the knob (9.4) must point to the center position (9.5). The width of the dividing line is then set with an Allen key situated in a holder on the back of the stand.

Caution!

Remember to put the Allen key back in the holder for safekeeping after you have set the dividing line.

It is especially important that the Allen key is not left in the adjustment opening when the knob (9.4) for positioning the dividing line is rotated.

8.11 Magnification calibration

There is no need to match the right and left magnification as is necessary, for example, for zoom-based optics, nor to calibrate the optics with a stage micrometer, as the maximum variance is 0.4%.

9. Accessories

A specially developed work table and correctly sized work chair allow ergonomically optimized work.

In addition to various optical accessories, such as filters and stage scales, a number of different slide holders can be used.

9.1 Motorized work table

The work table with motorized height adjustment is designed for optimum ergonomic conditions when working at the comparison microscope. Its table surface measures 120 cm x 56 cm and it is mounted symmetrically on the height-adjustable column.

Caution!

When setting up the comparison microscope, make sure that the stand is positioned exactly over the adjustable-height column.

Fig. 30 Motorized work table

- 1 Buttons for height adjustment
- 2 Adjustable-height feet

After the power supply is connected, the work table can be continuously adjusted within a range of approx. 62 to 92 cm over floor height. To do so, use the buttons (30.1).

The lifting speed is approx. 12 mm/sec. The lifting capacity is approx. 2000 N.

If the floor on which the work table is installed is not completely level, differences in level can be compensated for by screwing the feet (30.2) in or out

9.2 Stackable stages

The macroscopes are equipped with manual cross-stages.

Stage surface: 140 mm x 202 mm
Travel range of the stage: 50 mm x 50 mm

On the stage surface, there are threaded bore holes as receptacles for the adjustable holder (universal holder), rotating stages \varnothing 118 mm, large sample stages and rotation device for tool-mark examinations.


Rotating stages \varnothing 118 mm

Rotating stages, which can be stacked on the basic equipment cross-stages (Fig. 32), are fastened using two Allen screws. They are suitable for incident light and transmitted light examinations. The rotating stages each have one opening in the center, which is sealed with a removable glass plate (32.3). The threaded bore holes on the stage surfaces are provided for screwing on the adjustable holder. Additionally, there are also holes for inserting the stage clips. The rotating stage surface can be locked into place using a fastening screw (32.2).

In the bottom parts of the rotating stages, there are receptacles for inserting the polarizers (32.1) with λ -plates (red = OK). These can be rotated 90° .

Fig. 32 Rotating stage \varnothing 118 mm

- 1 Filter slot for polarizer
- 2 Fastening screw
- 3 Glass plate

Large sample stages

The stackable large sample stages can be screwed directly onto the cross-stages. They are primarily used for holding documents. Documents can be securely held in place and stretched using the magnetic clamps provided (33.1).

Of course, the large sample stages can also be used to hold other large-surface carrier material, such as materials bearing toolmarks.

Inclining rotating stages \varnothing 75 mm

The inclining rotating stages are likewise screwed directly onto the cross-stages. After releasing the lever (34.1), you can rotate them horizontally and incline them up to 45° using the built-in ball-and-socket joint (34.2). After positioning the comparison specimens, be sure to close the clamping system using the lever.

The surface of the stages (which are 75 mm in diameter) has a ridged surface (34.3) for fastening specimens using plasticine. Four bore holes are provided, into which stage clips (34.4) can be inserted.

In the center of the stage plate, there is a conical bore hole (34.5), into which all specimen receptacles can be inserted (in a similar fashion as the adjustable holder)

Fig. 33 Large sample stage (assembly)

1 Magnetic clamps



Fig. 34 Inclining rotating stage \varnothing 75 mm

1 Clamping lever, 2 Ball-and-socket joint, 3 Ridged surface
4 Bore hole for stage clips, 5 Conical bore hole

9.3 Slide holder

Adjustable holder with integrated rotating stage

Fasten the adjustable holders (universal holders) (35.1) to the cross-stage using the two screws (35.4).

After the clamping screw (35.3) is released, the universal holder can be moved sideways. This is especially useful for long, horizontally clamped pieces of evidence.

You can adjust the adapter sleeve (35.6) in the segment guide to an angle of up to 90° after releasing the clamping screw (35.2).

The adapter sleeve itself has a rotation device with gradation (35.7) and can be fixed into place using the clamping screw. The receptacles for pieces of evidence listed below can be inserted into the conical bore of the adapter sleeve (35.6).

After the clamping screw (35.5) is unscrewed, the entire unit can be rotated around its own axis.

Fig. 35 Adjustable holder

- 1 Adjustable holder
- 2 Locking device: angle adjustment of the adapter sleeve
- 3 Locking device: lateral adjustment of the holder
- 4 Fastening onto stage
- 5 Locking device: rotation of the holder
- 6 Adapter sleeve
- 7 Rotation device with gradation and clamping screw

Bullet Holder

To assemble the bullet holder (36.1), slide the adapter sleeve in the segment guide downwards to the "0" marking. Guide the bullet holder vertically to the adapter sleeve and fasten it with the knurled screw (36.2). Adjust the adapter sleeve horizontally again (slide upwards to the "90" marking) and press the rubber-coated pressure plate (37.3) into the conical bore of the adapter sleeve. After unscrewing the clamping screw (37.2), you can move the centering insert (37.3) horizontally in the bullet holder, thus adjusting it to the length of the projectiles currently being compared.

Caution!

Before removing the unit, be sure there is sufficient distance between it and the objective.

Fig. 36 Adjustable holder with bullet holder (assembly)

- 1 Bullet holder
- 2 Fastening screw



Fig. 37 Adjustable holder with bullet holder

- 1 Rubber-coated pressure plate
- 2 Clamping screw
- 3 Centering insert with spring system
- 4 Bullet

Paste-on Cylinder for Bullets

Using adhesive – or using piceine – you can stick (for example) deformed bullets that fit neither into the bullet receptacles nor into the receptacles for deformed pieces of ammunition onto the paste-on cylinder (38.1). There are paste-on cylinders with diameters of 5, 6, 8 and 10 mm.

Receptacles for deformed bullets

Receptacles for deformed bullets (39.1) are available for holding deformed bullets, which cannot be held by the receptacles described above. Their cavities sufficiently hold such specimens in place.

Fig. 38 Adjustable holder with paste-on cylinder for bullets

- 1 Paste-on cylinder
- 2 Bullet

Fig. 39 Adjustable holder with receptacle for deformed bullets

- 1 Receptacle for deformed bullets

Plug Receptacles for Bullet Casings

For fastening undamaged bullet casings of various calibers, plug receptacles (40.1) are provided. Their plastic holding devices are perfectly suited for holding bullet casings (40.2). The various calibers are marked on the rear sides of the plug receptacles. Worn-out holding devices can be replaced. Corresponding replacement parts are included in the delivery.

For deformed bullet casings, the brush holder set (40a) is available.

Plug receptacles for air rifle ammunition

The stainless steel pyramids (41.1) are intended for attaching air rifle ammunition.

Fig. 40a Brush holder set for bullet casings
(deformed bullet casings are also accommodated)

Fig. 40 Adjustable holder with plug receptacle for bullet casings

- 1 Plug receptacle
- 2 Bullet casing



Fig. 41 Adjustable holder with plug receptacles for air rifle ammunition

- 1 Plug receptacle

9. Accessories

Shotgun bullet casing holder

The shotgun bullet casing holder (42.1) has three spring rods rounded at the edges (42.2), to which shotgun bullet casings (43.1) of various calibers can be attached. The conical base (43.2) is provided for centering the bullet casings.

Fig. 42 Adjustable holder with shotgun bullet casing holder

- 1 Shotgun bullet casing holder
- 2 Spring rods



Fig. 43 Adjustable holder with shotgun bullet casing holder

- 1 Shotgun bullet casing
- 2 Conical base

Wire holder

In the wire holder (44.1), round parts between 0.6 mm and 8 mm in diameter can be clamped in order to (for example) examine cutting marks.

Receptacles for locking cylinders

Using an underlay of plasticine, these receptacles (45.1) can be used to securely fasten lock cylinders of safety locks bearing toolmarks.

Fig. 44 Adjustable holder with wire holder

- 1 Wire holder

Fig. 45 Adjustable holder with receptacle for locking cylinders

- 1 Receptacle for locking cylinder

Stages ø 20 mm

You can fasten pieces of evidence of various dimensions to the ridged surface of the stages with ø 20 mm (46.1). As the adhesive, you can use plasticine or wax.

Stages ø 60 mm

Specimens with a circumference of up to 55 mm and a thickness of up to 15 mm can be fastened with clamps (47.2) to the stages ø 60 mm (47.1). To do so, release the screws (47.3), raise the clamps, attach the specimens and press the clamps back in before tightening the screws. The additional threading (47.4) on the border of the stages ø 60 mm is for arranging the clamps asymmetrically. If you remove these, the stages ø 60 mm can also be used as large paste-on stages. You can fasten pieces of evidence of various dimensions to the ridged surface of the stages using plasticine.

Fig. 46 Adjustable holder with stage ø 20 mm

- 1 Stage ø 20 mm
- 2 Joint holder



Fig. 47 Adjustable holder with stage ø 60 mm

- 1 Stage ø 60 mm
- 2 Clamp
- 3 Fastening screw
- 4 Thread

Stages ø 25 mm with ball-and-socket joint

The stages ø 25 mm with ball-and-socket joint (48.1) are inserted into the adapter sleeve of the adjustable holder.

After releasing the clamping screw (48.2) on the ball-and-socket joint, the stage surface, which has a ridged surface, can be inclined arbitrarily. After aligning the specimen, retighten the clamping screw (48.2).

Joint holder

If, for example, you want to compare the casing sleeves of bullet casings, the profiles of the strike pin impressions on the percussion caps or the impressions of the extractor one after the other in one clamping, it is necessary to use joint holders (49.1). Its joint is fastened with a knurled screw (49.2).

The joint holders are provided for use as connecting pieces for the above-mentioned receptacles, with the exception of the bullet holder. In addition to the adjustment options provided by the adjustable holders, they offer an additional option for rotating and inclining the specimens. The incline of the specimens can be locked in place using the knurled screw (49.2)


Fig. 48 Adjustable holder with stage ø 25 mm with ball-and-socket joint

- 1 Stage ø 25 mm with ball-and-socket joint
- 2 Clamping screw

Fig. 49 Adjustable holder with joint holder

- 1 Joint holder
- 2 Clamping screw

9.4 Universal rotating mechanism

The rotating mechanism (order No. 581 088) (Fig. 52) consists of a stage plate, which can be mounted directly onto the cross-stage, and an articulated illumination arm, which allows you to rotate the illumination around the specimen.


Either the inclining rotating stage \varnothing 75 mm or the universal slide holder can be installed on the rotating plate..

Fig. 52 Universal rotating mechanism

Fig. 53a Light panel

1 On/off switch


9.5 Light panels

For homogeneous illumination of larger specimens, you can use the light panels (Fig. 53a, 53b). The light panel (Fig. 53a) comes with an articulated arm which can be attached to the universal object holder or the stage bracket.

Alternatively, the light panel shown in Fig. 53b can be used with the cold light illuminator and a fiber optic guide.

Fig. 53b Light panel

9.6 Filter for color contrast

For the differentiation of color deviations of two specimens being compared, the two filters CRA (red) and CGA (green) in sliders (Fig. 54) are provided. They are inserted into the filter receptacles in the macroscope holders (55.1) above the objectives.

For adjusting the superimposed image with color contrast, first set the composite image method **MIX**.

If the specimens of comparison overlap in their essential structures, these only appear in their original color where matching details are present.

Depending on the assignments of the filters to the two beam paths, the deviations are depicted in either red or green.

Fig. 54 Color filter

Fig. 55 Objective turret with objectives

1 Filter receptacle

9.7 Further illumination options

As well as the LED spot light (→ p. 14) and the light panels (→ p. 31) there are many other illumination options.

LED illumination

- LED power spot with one LED (3.5 W)
- LED spot with 19 LEDs
- LED ring illuminator (Fig. 57) with partial switching of the LED segments for close-up vertical illumination

To set the intensity of the illumination, refer to the separate “LED Controller Coolcontrol” manual provided.

UV illumination

The UV ring illuminator is available for UV illumination.

Cold light illumination

The cold light illuminator, e.g. KL2500 LCD (Fig. 56) can be installed next to the stand or on a mounting plate that is screwed onto the base of the Comparison Microscope.

The illumination is provided via a two-arm flexible fiber optic guide.

The cold light illuminator is controlled by the rotary knobs (56.2).

Fig. 56 Cold light illuminator, e.g. KL2500 LCD

- 1 Connection for light guide
- 2 Controls
- 3 Mounting plate


Fig. 57 LED ring illuminator

Fig. 58 UV ring illuminator

9.8 Stage scales and calibration standard

During assembly of the comparison microscope, the macro-objectives have been chosen so that they match with regard to their optical properties and permissible tolerances. For identical objectives, the image scales generated thus fall within the permissible tolerances.

The stage scale (stage micrometer) 11 519 963 or the calibration standard 11 581 081 can be used to verify this. To do so, adjust the section image for identical specimen magnifications and check that the image scales on both sides match exactly.

Fig. 59 Stage micrometer

10. Technical Data

Stand

- Stable and torsion-free cast-iron stand with motorized height adjustment for the carrier of the comparison bridge (stroke 255 mm)

Focusing

- Two focusing drives
- Stroke: 25 mm

Stages

- cross-stages
- Stage surface: 140 mm x 202 mm
- Travel range: 50 mm x 50 mm
- Ergonomic operating knobs for cross and synchronized displacement
- Bore hole receptacles for the adjustable holder, the rotating stages, the large sample stages or the bullet holder

Stackable rotating stages ø 118 mm

- Openings ø 50 mm
- Removable glass stage plates
- Clamping device for the stage rotation
- Stage carrier with receptacle for rotating polarizer

Inclining stages ø 75 mm

- Can be locked into place
- Inclines in all directions up to 45°
- Surface with ridged surface

Large sample stages

- Stage surface: 210 mm x 300 mm
- Metal plates for thin specimens such as documents
- Can be stacked on the cross-stage

Comparison bridge with ergonomic tube

- For superimposed image or split image comparisons with variable dividing line
- Variable width of dividing line
- Combination of split image and superimposed image possible
- Color variation of deviation in traces possible in superimposed image observations
- Maximum variance between right and left magnification: 0.4%
- Receptacles for macro-objectives 0.4x, 1x, 2x, 4x, 8x with adjustable iris diaphragms
- Rotating fasteners for the holders for oblique incident illumination
- Slots for inserting filter slides
- Distance between the optical axes: 450 mm

Tube factor

1x, 1.5x, 2x with magnification changer

Field number

FOV 22

Image orientation

Portrait or landscape

10. Technical Data

Objectives

Macro-objectives 0.4x, 1x, 2x, 4x, 8x

Eyepieces

HC PLAN S 10x/22

HC PLAN S 10x/20

HC PLAN S 10x/18

Magnifications and object fields

Magnification changer in position	Objectives	Free working distance	SFZ with eyepiece 10x	Total magnification with eyepiece 10x
	0.33:1	100 mm	166.0 mm	1.32x
	0.4x	60 mm	55.0 mm	4x
1.5x	0.4x	60 mm	36.6 mm	6x
2x	0.4x	60 mm	27.5 mm	8x
	1x	60 mm	22.0 mm	10x
1.5x	1x	60 mm	14.6 mm	15x
	2x	60 mm	11.0 mm	20x
1.5x	2x	60 mm	7.3 mm	30x
	4x	60 mm	5.5 mm	40x
1.5x	4x	60 mm	3.6 mm	60x
	8x	48 mm	2.75 mm	80x
1.5x	8x	48 mm	1.8 mm	120x
2x	8x	48 mm	1.37 mm	160x

Motorized work table

- Height adjustment range: 619 mm + 300 mm movement
- Lifting capacity: 2000 N
- Lifting speed: approx. 12 mm/s
- Load capacity: 200 kg
- Control system: up and down switches and foot switch
- Stage surface: 1200 mm x 560 mm
- Universal power supply

General specifications

- For indoor use only.
- | | |
|-----------------------|------------------------|
| Supply voltage: | 100-240 VAC |
| Frequency: | 50-60 Hz |
| Power input: | I _{max} 0.5 A |
| Fuse (power supply): | 2 x 2.5 A H 250 V |
| Ambient temperature: | 15-35°C |
| Relative humidity: | max. 80% up to 30°C |
| Overvoltage category: | II |
| Pollution degree: | 2 |

Dimensions without camera

- Height: 845 mm
(maximum; Z column extended)
- Width: 700 mm
(x/y stages extended to maximum width)
- Depth: 530 mm
(incl. front control knobs)
- Viewing height for ballistic objects in the universal holder: 590 mm – 510 mm

Weights

- | | |
|-------------------------------|-------|
| Basic equipment: | 45 kg |
| Comparison bridge: | 15 kg |
| Macroscopic stand with table: | 32 kg |
| Mot. work table: | 50 kg |

11. Care of the Instrument

Caution!

Unplug the power supply before performing cleaning and maintenance work!
Protect electrical components from moisture!

Microscopes in warm and warm-damp climates require special care in order to prevent the build up of fungus.

The instrument should be cleaned after each use, and the optics must be kept extremely clean.

11.1 Dust cover

Note:

To protect against dust, cover the instrument and accessories with the dust cover after each use.

11.2 Cleaning

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with little soapy water, benzine or ethyl alcohol.

For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Acetone, xylene or nitro-containing thinners can harm the instrument and thus must not be used.

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure when doing so that coated or plastic surfaces do not become matted or etched.

Cleaning the Stage

Remove light-colored spots on the stage by rubbing with paraffin oil or acid-free Vaseline.

Cleaning Glass Surfaces

Remove dust on glass surfaces with a fine, dry and grease-free hair brush, by blowing with a bellows or vacuum suction.

Carefully remove stubborn dirt on glass surfaces with a clean cloth moistened with distilled water. If the dirt still cannot be removed, use pure ethyl alcohol or benzine.

Cleaning Objectives

Caution!

The objective may not be unscrewed during cleaning. If damage appears on inner surfaces, the objectives must be sent to your Leica subsidiary for repair. We also advise against cleaning the inside surfaces of the eyepieces.

The front lenses of objectives are cleaned as described under "Cleaning Glass Surfaces". The upper lens is cleaned by being blown off with a pneumatic pump.

Removing Immersion Oil

Caution!

Follow safety instructions for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

11.3 Handling acids and alkalines

For examinations using acids or other aggressive chemicals, particular caution must be taken.

! Caution:

Be absolutely certain to prevent the optics and mechanical parts from coming into direct contact with these chemicals.

12. Essential Wear and Spare Parts

Order No. Material No.	Name	Used for
<u>Fuses:</u> 11362150001201	2 x 2,5 A H 250 V	Fuse (power supply)

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